



ZAMMAR[®]



A Revolution in Temperature Treatment



Countless advantages no contra-indications

- Continuous programmable temperature (cold and heat).
- Highly effective circulatory and muscle relaxant properties.
- Extremely easy to use both during competitions and in the stable.





-
- Avoids the administration of anti-inflammatory and other drugs with their undesired side effects.
 - Safe and efficacious treatment, highly innovative compared to empirical methods (irrigation, ice or glyco-gels).
 - Our technical department is at your disposal to design wraps suitable for special needs.
-



NEW

STABLE



Particularly suitable for stables

“STABLE” is a new product, developed on the basis of experience in the field, which offers:

- Increased refrigerating capacity
- **Possibility of housing all accessories necessary for treatments**
- Easy transport (on wheels)
- Strength (steel plate, treated and painted)
- Anti-dust filters to prevent obstructions in condenser

ACCESSORIES



Specially adapted folding trolley, compact and easy to handling. Protection curtain for unit

Support to be applied to the box wall, so that the unit is kept in a safer place and is less exposed to dust (decreased possibility of obstructions).



Transport bag for pipes and wraps

Connecting hoses equipped with a patented safety system which, in case of escape, releases the animal without hurting it or damaging the unit.



Models











**PATENT
PENDING**

A proven product



	STANDARD		POWERFULL		STABLE
Model	MS 365		MS 385		MS 599
Version	AA	BA	AA	BA	AA
Voltage	230/1/50 Hz	230/1/50 Hz	230/1/50 Hz	230/1/50 Hz	230/1/50 Hz
Operating mode	●● cold/hot	● cold	●● cold/hot	● cold	●● cold/hot
Temperature range for set point cold mode	0 / +10	0 / +10	-3 / +10	-3 / +10	-5 / +10
Temperature range for set point hot mode	+30 / +40	-	+30 / +40	-	+30 / +45
Work / pause times	up to 60 min. MAX		up to 60 min. MAX		up to 60 min. MAX
Massage	Adjustable		Adjustable		Adjustable

Combination of Thermic wraps

 Neck	X	X	X
 Shoulder	X	X	X
 Tendon	XX	XXXX	XXXX
 Knee • Tendon • Fetlock (anatomic)	XX	XXXX	XXXX
 Knee • Tendon • Fetlock Hoof (anatomic)	-	XX	XX
 Hoof	XX	XXXX	XXXX
 Stifle	X	X	X
 Hock	XX	XX	XX
 Back	X	X	X
 Dorsal cover sheet	-	X	X

X=Nr Wraps

Above data refer to an outside temperature of +30 °C

Descriptions, technical data and pictures are subject to modifications. For technical and commercial reasons, ZAMAR reserves the right to make changes at any time without notice, as required by our product development policy.



www.zamar.it

Strada Naz. della Cisa Km162, 34/A
46029 Suzzara (MN) Italy
Tel. 0039 0376 522158
Fax 0039 0376 509454
E-mail: info@zamar.it

Authorized dealer

Ludger Beerbaum has adopted
the ZAMAR device



Sport...and more

ZAMAR is typically used in competitions, but its benefits can also be experienced in daily use. All riding disciplines can take advantage of this advanced method for the treatment of tendinous and muscular problems both before and after a competition as well as in daily treatment.

Strongly recommended by vets.

The ZAMAR device has been tested and used by vet clinics with top-level experience, which work in direct contact with the national and international competitive environment.

The results are highly satisfactory as for immediate therapeutic response, simplicity of use and versatility of the unit.


Countless pathological conditions due to trauma or extreme work of the muscle-skeletal apparatus can be found in every sport.

The ZAMAR device has been tested in many ways and at significant times in all riding sports.

Here below we supply a list of pathologies suffered by sports horses after particularly demanding competitions, which were treated with the ZAMAR machine:

- **muscle/tendon contractures and injuries,**
- **inflammation and swelling due to stress and strain,**
- **pain and heaviness due to muscle strain,**
- **lactic acid stagnation,**
- **haematoma and extravasation due to contusion, trauma and tissue damage.**

The results were excellent. Great appreciation has been expressed for the possibility of treating the injured/inflamed parts with personalised therapies: cold and heat treatment, alternated cold and hot cycles and massage.



A fundamental resource: steady temperature

Conventional treatments (ice packs or glycol bags) do not ensure a steady temperature during the treatment. **ZAMAR has developed a new method which represents a safe and definitive solution to the problem of cold/heat therapy.**